



ALBUMIN (BCG Dye Method-End Point)

Intended Use

The reagents are used for the quantitative determination of Albumin in serum or plasma. For in-vitro diagnostic use only.

Introduction

Albumin plays important physiological role in the body. In circulation of blood, 80% oncotic pressure is maintained by it. Albumin assists in the transportation of bilirubin, calcium, fatty acids, bile acids, and drugs. It sequesters the toxins entering in to the body fluid and delivers them to the liver for detoxification. Testing albumin levels in serum helps to assess a person's nutritional status and risk for malnutrition. The test helps in determining if a patient has liver disease or kidney disease, or if not enough protein is being absorbed by the body.

Method

Bromocresol green method (Dye Method- End point)

Principle

Albumin forms blue green colored complex with Bromocresol green in acidic medium. The color formed is directly proportional to the concentration of albumin present in the sample, intensity of color is measured at 620nm.

Reagent Composition

Reagent 1: BCG Reagent

| | |
|-------------------|-------------|
| Bromocresol green | 0.15 mmol/l |
| Succinate buffer | 75 mmol/l |

Reagent 2:

| | |
|------------------|--------|
| Albumin Standard | 3 g/dl |
|------------------|--------|

Precautions

Following precaution should be taken:

- Avoid ingestion, do not pipette by mouth.
- Avoid contact with skin and eyes. If spilled, thoroughly wash affected area with water.
- Flush with plenty of water while disposing.

Reagent Storage and Stability

Unopened Reagent 1 is stable till expiry mentioned on the label when stored at room temperature.

Standard Reagent R2 is stable till expiry mentioned on the label when stored at 2-8°C.

Note: On request, Reagent 3 (Albumin: 4 g/dl) & Reagent 4 (Albumin: 6 g/dl) can be provided for linearity check with Reagent 2 (Albumin: 3 g/dl-Standard).

Reagent Preparation

Reagents are ready for use.

Reagent Deterioration

Reagents should be clear solutions. Turbidity and/or precipitation may be because of deterioration.

Sample Collection and Storage

Serum or plasma (EDTA, citrate or heparin) collected as per the standard procedures.

Albumin in serum is stable for 3 days at 2-8°C

General Assay Parameters

| | |
|----------------------------------|-----------|
| Mode | End Point |
| Wavelength (nm) | 620 |
| Wavelength Range Usable (nm) | 570-650 |
| Sample Volume (µl) | 10 |
| Reagent R1 (µl) | 1000 |
| Incubation Time (min.) | 5 |
| Incubation Temp. (°C) | RT |
| Normal Low (g/dl) | 3.2 |
| Normal High (g/dl) | 5.0 |
| Linearity (g/dl) | Upto 6 |
| Concentration of Standard (g/dl) | 3 |
| Blank with | Reagent |
| Units | g/dl |

Procedure

One reagent blank and one standard are sufficient for each assay series.

Pipette into test tubes:

| | Blank | Standard | Sample |
|-------------|--------|----------|--------|
| Reagent 1 | 1000µl | 1000µl | 1000µl |
| Dist. Water | 10µl | - | - |
| Reagent 2 | - | 10µl | - |
| Sample | - | - | 10µl |

Mix well & incubate for 5 min at room temperature.

Measure the absorbance of standard (A std) and sample (A sample) against reagent blank at 620nm.

Calculation

Albumin concentration in the sample can be calculated using the following formula:

$$\text{Albumin (g/dl)} = \frac{\text{Absorbance of Sample}}{\text{Absorbance of Standard}} \times \text{Conc. of Std.}$$

Example: If the absorbance of sample is 0.200 and the absorbance of standard is 0.180. The calculation shall be:

$$\frac{0.200}{0.180} \times 3.0 = 3.33 \text{ g/dl}$$

If the albumin concentration exceeds 6 gm/dl, dilute the sample with normal saline and repeat the assay. Get the results by multiplying by dilution factor.

Expected values

Neonates: 3.8 – 4.2 gm/dl
 Adult : 3.8 – 4.4 gm/dl

Note:

Expected range varies from population to population; therefore, each laboratory should establish its own normal range.

Limitations

1. Lipemic samples may give false elevated results. Lipemic samples should be reported with serum blank as below:
 - a) Add 10µl of serum sample to 1000µl of DI water and read absorbance at 620 nm.
 - b) Subtract the absorbance obtained as above from the absorbance of sample. Use this corrected absorbance for calculation.
2. The reagent and sample volumes can be altered proportionately so that the sample:reagent ratio remains same.
3. Calibration using aqueous standard provided in the kit may cause a matrix effect in some analyzers. In such a case it is recommended to use a serum based calibrator.

Quality Control

The patient results obtained for each batch can be validated by using normal and abnormal control sera with assayed values for albumin.

Performance

Linearity limit: 6 g/dl

Precision

Within run

| Control | Control 1 | Control 2 |
|----------------|-----------|-----------|
| No. of samples | 20 | 20 |
| Mean (g/dl) | 3.83 | 4.25 |
| S.D. | 0.05 | 0.04 |
| C.V. % | 1.29 | 0.84 |

Between run

| Control | Control 1 | Control 2 |
|----------------|-----------|-----------|
| No. of samples | 60 | 60 |
| Mean (g/dl) | 3.71 | 4.23 |
| S.D. | 0.07 | 0.06 |
| C.V. % | 1.75 | 1.36 |

References

- 1 Dumas B. T. , Arends R. L. , Pinto P. C. in standard Methods of Clinical Chemistry (1972) Vol. 7, p. 175-189, Academic Press Chicago.
2. Tietz N. W. (Ed.), Textbook of Clinical Chemistry, W. B. Saunders (1986) p. 589.
3. Young DS. Effects of drugs on CLinical Laboratory Tests, AACC Press, Washington DC, 5th edn. 2000.

Pack Presentation

| Product Code/ Catalogue No. | Pack Size* | Reagent 1 | Reagent 2 |
|--------------------------------|---------------|-----------|-----------|
| KGALB101.1.1 | 2x50ml | 2x50ml | 1x2ml |
| KGALB101.1.2 | 5x50ml | 5x50ml | 1x3ml |











* Pack size may vary on market demand.

Revision No: (Ver: KGALB101.1/1)

Date of Issue:1st April 2010

Symbols

Following symbols are used in the labeling of KEE GAD kits:

| | | | |
|---|----------------------|---|---------------------|
|  | Catalogue No. |  | Batch No. |
|  | CE Mark |  | Read instructions |
|  | In Vitro Diagnostics |  | Storage temperature |
|  | Expiry Date |  | Content |
|  | Product Name |  | Manufactured By |



Manufactured by:
KEE GAD Biogen Pvt. Ltd.
 A-8, Third Floor, Naraina Industrial Area,
 Phase-II, New Delhi-110028 (India)